

SINOHERO SW100 is a reliable, smart and high performance device for warming infusions, blood and transfusions to prevent any condition of hypothermia. Device temperature can be set freely in increments of 0.1°C from 35°C to 42°C. The integrated universal mounting clamp can be placed on an infusion stand, IV pole as well as a normal rail.



Revolutionary Two Warming Mechanism

- IV set or feeding set can be warmed up 60cm in heating plate, much longer than other ordinary warmer with heating plate.
- The detachable tube warmer ensures warmed fluid up to the patient without heat loss, much more efficient.



Features

- Compatible for any standard I.V Set
- Easy and quick to set up temperature range 35~42°C
- Real-time display temperature
- · No additional disposable required
- · High optimized warming circuit
- Minimised Energy Consumption
- · Easy to clean



- Model: SW100
- · Device Name: Blood, Infusion, Feeding Warmer
- Compatible tube: Transfusion set, infusion set or feeding set with diameter 3.5-5mm
- Adapter Output: DC 15V, 2.0V
- Protection class: I
- · Heating board Degree of protection: BF
- Main engine Protection against ingress of water: IPX2
- Control temperature: 35°C 42°C adjustable in 0.1°C steps
- Temperature accuracy: ±0.5°C (heating plate), ±1°C (warming tube)
- Operational mode: Continuous operation
- Dimensions: 140x75x35mm
- · Weight: 210g
- Service life: 5 years at leastTube warmer Length: 800mm
- · Compact and user-friendly design for easy operation
- · Long warmer tube support flexible performance in multiple situations









What and why is tube warmer SW-WT?

SW-WT is a reusable silicone heated tube which connect and control by SW100, the temperature will show at the SW-100 display. This tube will highly improve the efficiency of heating, and optimisation between warmer and patient at low flow rates, so it is highly recommend option!

Specification subject to be changed without prior notice



Sino-Hero (Shenzhen) Bio-Medical Electronics Co., Ltd

Add: A area, Second Floor, First Building, Tongkangfu Industrial Park, Yingrenshi Community, Shiyan Sub-district, Bao'an District, 518108 Shenzhen, PEOPLE'S REPUBLIC OF CHINA Tel: +86-755-26408541 | 26404878 Fax: +86-755-26420566

Email: sales@sinoheromed.com

Web:www.sinoheromed.com