
X-Ray Upgrade Solution





Veterinary Flat Panel DR System V-ray30

Shenzhen Landwind Medical Co., Ltd Floor 8, Block A4, Nanshan iPack, No.1001 Xueyuan Road, Nanshan District, Shenzhen China. 518071 Tel: +86 755 8652 2242

Fax: +86 755 8393 5015 Email: Sales@landwindmedical.com Web: www.landwindmedical.com

DISTRIBUTOR:	
2019.06	V1.0

Specifications and appearance are subject to change without prior notice.








- Automatically collect animal information when radiographic through DICOM work list
- Shorten the input time of patient information according to APR automatically inputting preset image processing parameters, ROI, Markers, LUT, etc.
- Rich image editing tools including contrast, inversion, flip, rotation, ROI and window adjustment.
- The APR parameters of dog, cat, rabbit, bird and lizard are preset.

Sample1



Limiting Resolution: 3.93lp/mm Effective Pixel Matrix: 3328 X 3328 (Pixels) Active Area: 422.7mm X 422.7mm A/D: 14 bits


Stationary X-Ray bed

Tabletop Size: 1410mmX605mmX785mm Column height: 1360mm X-ray tube assembly lifting distance: 580mm Lateral suspension move range: 425mm Longitudinal suspension move range: 240mm

Dual Platform Operation

Fixed LCD touch screen operation, support various animals, different weights and different body parts examination

Image Acquisition DR Workstation remote control



Auxiliary diagnostic measurement tools



Sample2



Patient date managing system



Sample3